

Alternating Air Pressure Mattress User Manual


1. Safety

- a. Please read this User Manual carefully to avoid risk of electrical shock.
- b. Disconnect the power plug from the power socket if this product is not being used.
- c. Do not operate or place Air Pump near damp or wet areas.
- d. Do not place any foreign object into any of the openings.

2. Introduction

- a. This Air Mattress has 2 independent air chambers that provide alternating pressure relief and massaging effect of air circulating underneath the user's body to move blood and lymph in their correct direction when the user is not active enough to move on his/her own.
- b. This product should be used by people of all ages during times of prolonged immobility and serious illnesses.

3. Operation Guide



4. Installation Instruction

- a. Open the package and inspect all parts to ensure they are in good condition.
- b. Unfold the Air Mattress and lay it on top of the mattress. Slip the Anti-slip Pad underneath the mattress.
- c. Place the Air Pump on the floor or hang it on the bed frame.
- d. Connect the Air Tubes between the Outlet of the Pump and the Inlet of the Air Mattress.
- e. Ensure the Air Tubes are not tangled or bent to avoid obstruction to air flow from the Pump to the Air Mattress.
- f. Connect the Power Plug to electric outlet and turn the Power Switch to ON.

- g. Do not place any object on the Air Mattress and allow the Air Pump to inflate the Air Mattress for at least 20 minutes or until its fully inflated.
- h. Once the Air Mattress is fully inflated, lay bed sheet over the Air Mattress and it is ready for use.
- i. Adjust the Pressure Adjustment Knob to obtain optimum pressure level to suit individual comfort.

5. Product Specification

a. Air Pump

- Input Voltage: 220V/50Hz
- Power: 10 Watts
- Alternating Pressure Cycle: every 6 minutes
- Maximum Pressure: 1.8 PSI
- Maximum Volume: 8 Liters / minute

b. Air Mattress

- Material: Vinyl
- Size: 190 x 90 cm
- Weight: 1.5 Kg
- Maximum Load: 100 Kg

